

Certificate of product line quality

ClipTip Pipette Tips product family

Thermo Fisher Scientific, Inc. certifies that the Thermo Scientific™ ClipTip™ Pipette Tips product family has been manufactured in accordance with established manufacturing guidelines and product specifications. These products are subject to rigorous quality inspections and conform to the highest internally established quality requirements as well as external standards as listed below. The following quality control criteria apply to the aforementioned product family:

Item	Description
Laboratory use, non-regulated	Based on the Intended Use and/or Indicates for Use, this product family is not considered to be a medical device.
Materials	Products are produced from virgin resins and are absent of recycled resins. All products are produced with materials such as polypropylene compound (Low Retention additive) and polyethylene (filter). Our material suppliers do not use or intentionally incorporate the following agents into the materials used for production: plasticizers (softeners/phthalates), or latex.
Production	Products are manufactured in a controlled environment.
Functional testing	Products undergo continuous quality controls regarding product function; this includes functional testing with the instrument (sample-based).
Animal components	Tallow-derived materials used in the product fulfil the requirements laid down in the Regulations 1069/2009/EC, and 142/2011/EC, and the "Note for Guidance EMA/410/01, and as amended.
Bisphenol A (BPA) free	BPA has not been used in the manufacturing process of the plastic in question
Adenosine triphosphate (ATP) free	Products are tested for the presence of ATP using a luciferase-based detection assay. Products with ATP levels below the experimental baseline value of 10-13 are certified ATP-free.

Available in Fit for Purpose



ISO 9001
Certified

ISO 14001
Certified

ISO 13485
Certified

Ravi Gupta
Vice President/General Manager,
Liquid Handling Solutions
Laboratory Plastics Essentials

Luis Perez
Director, Quality Assurance
Laboratory Products

Thermo Fisher Scientific, Inc. Quality Management System is certified in accordance with international Quality Management standards. All sites are certified to one or more Quality Management standards; refer to the [Site ISO Certifications](#) for copies of relevant ISO Certifications.

Changes to manufacturing procedures, packaging, and product specifications are documented and traceable.

All products are fully traceable by lot number. For a lot-specific certificate, please visit us at thermofisher.com.

For any additional information or inquires, Contact Us for Support ([Link](#)) and choose the appropriate option.

Learn more at thermofisher.com/fitforpurpose

thermo scientific